# Direct Cosmetic RESTORATION KIT

For tooth preparation & contouring of direct restorative material

Direct Cosmetic RESTORATION KIT is designed under the guidance of Dr. Sushil Koirala with the purpose to achieve efficiency and accuracy during tooth preparation and finishing of direct cosmetic restorations with consideration of the new criteria protocols for maximum preservation of tooth structure.

## Features

- Speedy, simple and accurate tooth preparation can be achieved with 8 kinds of specially selected Diamond Points FG.
- For a wide range of direct cosmetic restoration works.
- Diamond Points FG are hygienically organized in an autoclavable aluminum bur block for convenient sterilization and replacement.

## Description and Intended uses

### Cavity Preparation

Diamond Points FG # 340S : for Occlusal Cavity Preparation and Veneer Margin Preparation

Diamond Points FG # MI-2P : for Minimal Intervention, preparing small cavities on pits, lingual and interdental surfaces

### Veneer Preparation

Diamond Points FG # 121 : for Depth Cutting

Diamond Points FG # 102R : for Labial Reducing

Diamond Points FG # F216 : for Gross Finishing, Creating Grooves, Cavity Beveling and Occlusal trimming

### Fine Finishing


Diamond Points FG # SF104R : for Finishing Interdental, Cervical and Flat surfaces

Diamond Points FG # SF215 : for Finishing Interdental and Cervical surfaces

Diamond Points FG # SF219 : for Finishing Occlusal and Lingual surfaces

## Maximum recommended speed

- Less than 450,000 min<sup>-1</sup>: 340S
- Less than 300,000 min<sup>-1</sup>: MI-2P, 121, 102R, F216, SF104R, SF215
- Less than 160,000 min<sup>-1</sup>: SF219

## Note

- F : denotes Fine Grit
- SF : denotes SuperFine Grit

## Caution

- Operate each instrument within the maximum permissible speed only.
- Insert the shank completely into the chuck of the handpiece.
- Rotate each instrument prior to use to check safety and smooth operation.
- Operate at recommended pressure only.
- Do not use at an acute angle.
- Wear eye-protecting glasses when using such instrument.
- Use these products for indicated purpose only.
- Do not use instruments with any of the defects or damages including deformation, run-out, rust, surface flaw, distortion and taint.
- Use with water spray and continuous gentle pressure.
- For reuse of Diamond Points FG, remove all fouling with cleaning liquid or disinfectant and then sterilize by any of such methods: autoclave, chemiclave, dry-heat sterilization, or disinfect with chemicals, following the instructions for use by manufacturers of sterilization/disinfection chemicals. However, chlorine disinfectant must not be used to avoid erosion.

### MANUFACTURED BY

SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan

SHOFU INC. Singapore Branch 10 Science Park Road #03-12, The Alpha, Science Park II, Singapore 117684

SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA

SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8JU, UK

SHOFU DENTAL SUPPLIES (Shanghai) CO. LTD. No.645 Jiye Road, Sheshan Industry Park, Songjiang 201602, Shanghai, China

EU REPRESENTATIVE: SHOFU DENTAL GmbH Am Brühl 17, 40878 Ratingen, Germany

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
#### MANUFACTURED BY

SHOFU INC. 11 Kamitakamatsu-cho, Fukuine, Higashiyama-ku, Kyoto 605-0983, Japan  
SHOFU INC. Singapore Branch 10 Science Park Road #03-12, The Alpha, Science Park II, Singapore 117684  
SHOFU DENTAL CORPORATION 1225 Stone Drive, San Marcos, CA 92078-4059, USA  
SHOFU DENTAL PRODUCTS LTD. Duke's Factory, Chiddingstone Causeway, Tonbridge, Kent, TN11 8 JU, UK  
SHOFU DENTAL SUPPLIES (Shanghai) CO. LTD. No.645 Jiye Road, Sheshan Industry Park, Songjiang 201602, Shanghai, China  
EU REPRESENTATIVE: SHOFU DENTAL GmbH Am Brühl 17, 40878 Ratingen, Germany